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ENGLISH

ATTENTION

IMPORTANT INSTRUCTIONS. KEEP THEM FOR FUTURE USE.

- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliances by a person responsible for their safety.

- Keep plastic bags and other packaging elements out of the reach of children. They may be potentially dangerous.
- Do not plug in before making sure the voltage on the rating plate and that of your household are the same.
- Do not submerge or wet the appliance.
- Don't allow children play with the appliance. Keep the appliance out of the reach of children.
- () Do not use it wrinkled or folded.
- () Do not insert pins or sharp objects.
- () Not suitable for children under 3 years.
- Do not use if it is wet.
- Do not use it with disabled people, children or people not sensitive to heat.
- It must not be used with animals.
- Use for extended periods of time in its maximum position may cause skin burns.
- Do not sleep with the pillow turned on. If you are going to use it for a long period of time, the on switch must be placed in the position 1 (minimum heat).
- Always use it with the pillow slip in place.
- Do not place the knob (1) or the power cable underneath or on top of the pillow when it is turned on.
- Before storage, let it cool before folding. If you fold it, you must not leave objects on top of it, not even clothing.
- Examine it frequently to see whether there are any signs of damage or wear and tear. If such signs appear or if the apparatus has been misused, it must be checked by an Authorised Technician before further use.

- This apparatus is not intended for medical use in hospitals.
- It must not be used by children unless the handset has been pre-adjusted by a responsible person, or the child has at least been told how to control it properly and safely.
- When using for the first time, you may notice a slight odour due to the synthetic materials.
- Checking the pillow outside with your hands must not be taken into account as this may give a misleading result of the heat emitted.
- The cleaning and maintenance done by users, can not be done by children without any supervision.
- If the cable were to deteriorate it must be changed by an Authorised Service Centre.

INSTRUCTIONS OF USE

- Connect the appliance to the mains and, first of all, place the On switch (2) in its maximum position (no. 4). Keep it in this position for some minutes. Then, if you wish, you may turn it down.
- When turning the appliance on, the indicators (3) will turn green for a few seconds and then yellow. After a period of time the green and yellow lights will alternate continuously, which indicates connections and disconnections to keep the selected temperature.
- Notice the switch positions.

Position 1: Minimum heat.

Position 2-3: Medium heat.

Position 4: Maximum heat.

Position 0: Turned off.

- Stretch it well out on the desired part of your body, without folds or bends. Cover it with a towel, blanket or similar item. This will prevent loss of heat. It will be more effective when the pillow and its slip are in contact with your skin.

- The model PHC50 has been especially designed to relieve back and neck pains. The adhesive belt in the top allows you to adjust it in the neck and with the plastic band you can adjust it to the body. Nevertheless you can use this model in any other part of your body with the same therapeutic benefits.
- **WARNING:** If the indicator lights turns to red color, immediately disconnect the appliance from the mains. It must be checked by an Authorised Technician. This situation can be produced due to a not correct using, for example when it is folded or bent. When this happen, switch off the unit from the mains, stretch it completely and switch it on again. Once the red color indicator light has disappeared the appliance is correct but if the indicator light still shows the red color, the appliance must be checked.

AUTO OFF

- The pillow has a timer (AUTO OFF) that automatically disconnects it after 90 minutes have elapsed since starting to be used.
- If you wish to continue using it when such system has disconnected it, you must place the On switch (2) in the "0" position, wait a few seconds and then place the switch in the desired position.

MAINTENANCE AND CLEANING

- Before cleaning, make sure that it is disconnected from the mains and left to cool.
- The pillow slip may be washed in a washing-machine at no more than 40°C. Bleach must not be used nor must it be dried in a dryer. It must not be dry cleaned or ironed either.
- Do not wash the pillow or dry clean it. You can clean it with a slightly moistened sponge.
- You must avoid wetting the switch or handset when washing.
- When drying, the cable must be placed away from any water to prevent it from entering the handset.
- Do not use chemical or abrasive products, metallic scourers to clean the appliance.

ENVIRONMENT PROTECTION

- Dispose the appliance according to the local garbage regulations.
- Never throw it away to the rubbish. So you will help to the improvement of the environment.

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Resguardo para enviar a JATA por el S.A.T. como prueba de garantía en sus liquidaciones.
Electrodomésticos JATA, S.A.
Recibo S.A.T. N°.